

Table of contents

\\USER

<table border="1"><tr><td>PHANTOM</td></tr></table>	PHANTOM
PHANTOM	

<table border="1"><tr><td><table border="1"><tr><td>PHANTOM_TEST</td></tr></table></td></tr></table>	<table border="1"><tr><td>PHANTOM_TEST</td></tr></table>	PHANTOM_TEST
<table border="1"><tr><td>PHANTOM_TEST</td></tr></table>	PHANTOM_TEST	
PHANTOM_TEST		

<table border="1"><tr><td><table border="1"><tr><td>dcm2niix_issue870</td></tr></table></td></tr></table>	<table border="1"><tr><td>dcm2niix_issue870</td></tr></table>	dcm2niix_issue870
<table border="1"><tr><td>dcm2niix_issue870</td></tr></table>	dcm2niix_issue870	
dcm2niix_issue870		

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>loca_3plan_haste</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms2</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms5</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb2</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb5</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms2</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb2</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>loca_3plan_haste</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms2</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms5</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb2</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb5</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms2</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb2</td></tr></table></td></tr></table>	<table border="1"><tr><td>loca_3plan_haste</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms2</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms5</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb2</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb5</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms2</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb2</td></tr></table>	loca_3plan_haste	PRODUCT__ep2d_bold__p2_sms1	PRODUCT__ep2d_bold__p2_sms2	PRODUCT__ep2d_bold__p2_sms5	C2P__cmrr_mbep2d_bold__p2_mb1	C2P__cmrr_mbep2d_bold__p2_mb2	C2P__cmrr_mbep2d_bold__p2_mb5	PRODUCT__ep2d_diff__p2_sms1	PRODUCT__ep2d_diff__p2_sms2	C2P__cmrr_mbep2d_diff__p2_mb1	C2P__cmrr_mbep2d_diff__p2_mb2
<table border="1"><tr><td><table border="1"><tr><td>loca_3plan_haste</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms2</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms5</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb2</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb5</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms2</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb2</td></tr></table></td></tr></table>	<table border="1"><tr><td>loca_3plan_haste</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms2</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms5</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb2</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb5</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms2</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb2</td></tr></table>	loca_3plan_haste	PRODUCT__ep2d_bold__p2_sms1	PRODUCT__ep2d_bold__p2_sms2	PRODUCT__ep2d_bold__p2_sms5	C2P__cmrr_mbep2d_bold__p2_mb1	C2P__cmrr_mbep2d_bold__p2_mb2	C2P__cmrr_mbep2d_bold__p2_mb5	PRODUCT__ep2d_diff__p2_sms1	PRODUCT__ep2d_diff__p2_sms2	C2P__cmrr_mbep2d_diff__p2_mb1	C2P__cmrr_mbep2d_diff__p2_mb2	
<table border="1"><tr><td>loca_3plan_haste</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms2</td></tr><tr><td>PRODUCT__ep2d_bold__p2_sms5</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb2</td></tr><tr><td>C2P__cmrr_mbep2d_bold__p2_mb5</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms1</td></tr><tr><td>PRODUCT__ep2d_diff__p2_sms2</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb1</td></tr><tr><td>C2P__cmrr_mbep2d_diff__p2_mb2</td></tr></table>	loca_3plan_haste	PRODUCT__ep2d_bold__p2_sms1	PRODUCT__ep2d_bold__p2_sms2	PRODUCT__ep2d_bold__p2_sms5	C2P__cmrr_mbep2d_bold__p2_mb1	C2P__cmrr_mbep2d_bold__p2_mb2	C2P__cmrr_mbep2d_bold__p2_mb5	PRODUCT__ep2d_diff__p2_sms1	PRODUCT__ep2d_diff__p2_sms2	C2P__cmrr_mbep2d_diff__p2_mb1	C2P__cmrr_mbep2d_diff__p2_mb2		
loca_3plan_haste													
PRODUCT__ep2d_bold__p2_sms1													
PRODUCT__ep2d_bold__p2_sms2													
PRODUCT__ep2d_bold__p2_sms5													
C2P__cmrr_mbep2d_bold__p2_mb1													
C2P__cmrr_mbep2d_bold__p2_mb2													
C2P__cmrr_mbep2d_bold__p2_mb5													
PRODUCT__ep2d_diff__p2_sms1													
PRODUCT__ep2d_diff__p2_sms2													
C2P__cmrr_mbep2d_diff__p2_mb1													
C2P__cmrr_mbep2d_diff__p2_mb2													

\\USER\PHANTOM\PHANTOM_TEST\ldcm2niix_issue870\loca_3plan_haste

TA: 21 sec Coil Selection: Auto Voxel Size: 0.5×0.5×4.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	On
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	9
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	2
Slices	9
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	3
Slices	9
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	320 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	1000.0 ms
TE	72.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	1000.0 ms
TE	72.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle Mode	Constant
Flip Angle 1	140 deg

Contrast - Common

Flip Angle 2	90 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
Wrap-up Magn.	None
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	320 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
Base Resolution	320
Phase Resolution	100 %
Interpolation	2.00

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	TSE/Separate
Acceleration Factor PE	3
Reference Lines PE	24
Deep Resolve	Off
Phase Partial Fourier	5/8

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Off
Image Filter	Off

Geometry - Common

Slice Group	1
Slices	9
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	2
Slices	9
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	3
Slices	9

Geometry - Common

Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	320 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	1000.0 ms
Multi-Slice Mode	Single Shot
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	2
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	3
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Miscellaneous

Coil Focus	Flat
------------	------

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
---------	----------

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1000.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	320 mm
FOV Phase	100.0 %
Phase Resolution	100 %
Motion Correction	None

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	None
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing

Inline Composing	Off
------------------	-----

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	h
Dimension	2D
RF Pulse Type	Normal
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	488 Hz/Px
Echo Spacing	6.04 ms
Turbo Factor	320

Sequence - Part 2

Introduction	On
Motion Correction	None

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\USER\PHANTOM\PHANTOM_TEST\dcn2niix_issue870\PRODUCT__ep2d_bold__p2_sms1

TA: 14 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	33.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	33.00 ms
MTC	Off
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	7
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	2
Reference Lines PE	32
SMS Factor	1
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	180.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0

BOLD

Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	7
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	1042 Hz/Px
Echo Spacing	1.06 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
Ghost Reduction	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM_TEST\dcn2niix_issue870\PRODUCT__ep2d_bold__p2_sms2

TA: 14 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 4 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	36.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	36.00 ms
MTC	Off
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	7
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	2
Reference Lines PE	32
SMS Factor	2
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	180.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0

BOLD

Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	7
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	1042 Hz/Px
Echo Spacing	1.16 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
Ghost Reduction	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM_TEST\dcm2niix_issue870\PRODUCT__ep2d_bold__p2_sms5

TA: 14 sec Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc:: 10 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	39.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	39.00 ms
MTC	Off
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	7
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	2
Reference Lines PE	32
SMS Factor	5
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	180.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0

BOLD

Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	7
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	1042 Hz/Px
Echo Spacing	1.26 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
Ghost Reduction	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM_TEST\dcn2niix_issue870\C2P_cmrr_mbep2d_bold_p2_mb1

TA: 14 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	33.00 ms
Averages	1
Multi-band accel. factor	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	33.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	7
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Single-shot
Acceleration Factor PE	2
Reference Lines PE	32
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	180.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Multi-band accel. factor	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	7
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Dimension	2D
Excitation	Standard
Gradient Mode	Performance
Flow Compensation	None
Bandwidth	1000 Hz/Px
Echo Spacing	1.08 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
RF Spoiling	Off

Sequence - Special

Excite pulse duration	2560 us
Min. prep scans	0
Delay before PC scans	0 us
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Sinc exc. pulse BWTP	5.20
Physio recording	Off
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM_TEST\dcn2niix_issue870\C2P_cmrr_mbep2d_bold_p2_mb2

TA: 22 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	34.80 ms
Averages	1
Multi-band accel. factor	2
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	34.80 ms
MTC	Off
Magn. Preparation	None
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	7
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Single-shot
Acceleration Factor PE	2
Reference Lines PE	32
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	-180.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Multi-band accel. factor	2

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	7
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Dimension	2D
Excitation	Standard
Gradient Mode	Performance
Flow Compensation	None
Bandwidth	1000 Hz/Px
Echo Spacing	1.16 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
RF Spoiling	Off

Sequence - Special

Excite pulse duration	2560 us
Min. prep scans	0
Min. prep scans SB	0
Delay before PC scans	0 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
Opt. MB RF pulse BW	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Sinc exc. pulse BWTP	5.20
Physio recording	Off

Sequence - Special

Triggering scheme	Standard
-------------------	----------

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM_TEST\dcn2niix_issue870\C2P_cmrr_mbep2d_bold_p2_mb5

TA: 37 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	38.00 ms
Averages	1
Multi-band accel. factor	5
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	38.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	7
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Single-shot
Acceleration Factor PE	2
Reference Lines PE	32
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	5

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	-180.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Multi-band accel. factor	5

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	7
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Dimension	2D
Excitation	Standard
Gradient Mode	Performance
Flow Compensation	None
Bandwidth	1000 Hz/Px
Echo Spacing	1.26 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
RF Spoiling	Off

Sequence - Special

Excite pulse duration	3200 us
Min. prep scans	0
Min. prep scans SB	0
Delay before PC scans	0 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
Opt. MB RF pulse BW	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Sinc exc. pulse BWTP	5.20

Sequence - Special

Physio recording	Off
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM_TEST\dcn2niix_issue870\PRODUCT__ep2d_diff_p2_sms1

TA: 21 sec Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	80.00 ms
Concatenations	1
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	80.00 ms
MTC	Off
Magn. Preparation	None
Fat-Water Contrast	Fat Saturation
Fat Saturation	Weak
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %

Resolution - Common

Interpolation	Off
---------------	-----

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	2
Reference Lines PE	24
SMS Factor	1
Deep Resolve	Off
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	180.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Concatenations	1

Physio - PACE

Resp. Control	Off
Concatenations	1

Diff

Diffusion Mode	3-Scan Trace
Diff. Directions	3
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	2000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	Off
Trace Weighted Images	On
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Fast
Bandwidth	1000 Hz/Px
Echo Spacing	1.08 ms
Free Echo Spacing	Off
Optimization	None
EPI Factor	100

Sequence - Part 2

Introduction	Off
Phase Correction	Internal
Ghost Reduction	Off

Sequence - Assistant

SAR Assistant	Off
Optimization	None

\\USER\PHANTOM\PHANTOM_TEST\dcn2niix_issue870\PRODUCT__ep2d_diff_p2_sms2

TA: 18 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 4 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	80.00 ms
Concatenations	1
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	80.00 ms
MTC	Off
Magn. Preparation	None
Fat-Water Contrast	Fat Saturation
Fat Saturation	Weak
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %

Resolution - Common

Interpolation	Off
---------------	-----

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	2
Reference Lines PE	24
SMS Factor	2
Deep Resolve	Off
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	P >> A
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	180.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Concatenations	1

Physio - PACE

Resp. Control	Off
Concatenations	1

Diff

Diffusion Mode	3-Scan Trace
Diff. Directions	3
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	2000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	Off
Trace Weighted Images	On
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Fast
Bandwidth	1000 Hz/Px
Echo Spacing	1.18 ms
Free Echo Spacing	Off
Optimization	None
EPI Factor	100

Sequence - Part 2

Introduction	Off
Phase Correction	Internal
Ghost Reduction	Off

Sequence - Assistant

SAR Assistant	Off
Optimization	None

\\USER\PHANTOM\PHANTOM_TEST\dcmm2niix_issue870\C2P_cmrr_mbep2d_diff_p2_mb1

TA: 30 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	100.00 ms
Multi-band accel. factor	1
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	100.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	90 deg
Refocus flip angle	180 deg
Fat-Water Contrast	Fat Saturation
Grad. rev. fat suppr.	Enabled
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Single-shot
Acceleration Factor PE	2
Reference Lines PE	24
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Multi-band accel. factor	1

Physio - PACE

Resp. Control	Off
Multi-band accel. factor	1

Diff

Diffusion Mode	MDDW
Diff. Directions	6
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	2000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Dimension	2D
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Performance
Bandwidth	1000 Hz/Px
Echo Spacing	1.08 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
RF Spoiling	Off
Phase Correction	Internal

Sequence - Special

Min. prep scans	0
Delay before PC scans	0 us
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz

Sequence - Special

Physio recording	Off
------------------	-----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM_TEST\dcmm2niix_issue870\C2P_cmrr_mbep2d_diff_p2_mb2

TA: 51 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	100.00 ms
Multi-band accel. factor	2
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	100.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	90 deg
Refocus flip angle	180 deg
Fat-Water Contrast	Fat Saturation
Grad. rev. fat suppr.	Enabled
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Single-shot
Acceleration Factor PE	2
Reference Lines PE	24
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

Geometry - AutoAlign

Slice Group	1
Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A10.0 F35.0
L	0.0 mm
A	10.0 mm
F	35.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A10.0 F35.0 mm
Orientation	Transversal
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.240676 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Multi-band accel. factor	2

Physio - PACE

Resp. Control	Off
Multi-band accel. factor	2

Diff

Diffusion Mode	MDDW
Diff. Directions	6
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	2000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Dimension	2D
Excitation	Standard
Gradient Mode	Performance
Bandwidth	1000 Hz/Px
Echo Spacing	1.16 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
RF Spoiling	Off
Phase Correction	Internal

Sequence - Special

Excite pulse duration	2560 us
Refocus pulse duration	5120 us
Min. prep scans	0
Min. prep scans SB	0
Delay before PC scans	0 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
Opt. MB RF pulse BW	Off
Time-shifted MB RF	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
PF omits higher k-space	Off

Sequence - Special

Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Physio recording	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----