

Table of contents

\\USER

<table border="1"><tr><td>PHANTOM</td></tr></table>	PHANTOM
PHANTOM	

<table border="1"><tr><td><table border="1"><tr><td>PHANTOM</td></tr></table></td></tr></table>	<table border="1"><tr><td>PHANTOM</td></tr></table>	PHANTOM
<table border="1"><tr><td>PHANTOM</td></tr></table>	PHANTOM	
PHANTOM		

<table border="1"><tr><td><table border="1"><tr><td>check_MBfactor_dcm2niix_issue870</td></tr></table></td></tr></table>	<table border="1"><tr><td>check_MBfactor_dcm2niix_issue870</td></tr></table>	check_MBfactor_dcm2niix_issue870
<table border="1"><tr><td>check_MBfactor_dcm2niix_issue870</td></tr></table>	check_MBfactor_dcm2niix_issue870	
check_MBfactor_dcm2niix_issue870		

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>localizer</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>localizer</td></tr></table></td></tr></table>	<table border="1"><tr><td>localizer</td></tr></table>	localizer
<table border="1"><tr><td><table border="1"><tr><td>localizer</td></tr></table></td></tr></table>	<table border="1"><tr><td>localizer</td></tr></table>	localizer	
<table border="1"><tr><td>localizer</td></tr></table>	localizer		
localizer			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms1</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms1</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms1</td></tr></table>	PRODUCT__ep2d_bold__p3_sms1
<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms1</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms1</td></tr></table>	PRODUCT__ep2d_bold__p3_sms1	
<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms1</td></tr></table>	PRODUCT__ep2d_bold__p3_sms1		
PRODUCT__ep2d_bold__p3_sms1			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms2</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms2</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms2</td></tr></table>	PRODUCT__ep2d_bold__p3_sms2
<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms2</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms2</td></tr></table>	PRODUCT__ep2d_bold__p3_sms2	
<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms2</td></tr></table>	PRODUCT__ep2d_bold__p3_sms2		
PRODUCT__ep2d_bold__p3_sms2			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms5</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms5</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms5</td></tr></table>	PRODUCT__ep2d_bold__p3_sms5
<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms5</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms5</td></tr></table>	PRODUCT__ep2d_bold__p3_sms5	
<table border="1"><tr><td>PRODUCT__ep2d_bold__p3_sms5</td></tr></table>	PRODUCT__ep2d_bold__p3_sms5		
PRODUCT__ep2d_bold__p3_sms5			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb1</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb1</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb1</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb1
<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb1</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb1</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb1	
<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb1</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb1		
C2P__cmrr_mbep2d_bold__p3_mb1			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb2</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb2</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb2</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb2
<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb2</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb2</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb2	
<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb2</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb2		
C2P__cmrr_mbep2d_bold__p3_mb2			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb5</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb5</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb5</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb5
<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb5</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb5</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb5	
<table border="1"><tr><td>C2P__cmrr_mbep2d_bold__p3_mb5</td></tr></table>	C2P__cmrr_mbep2d_bold__p3_mb5		
C2P__cmrr_mbep2d_bold__p3_mb5			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms1</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms1</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms1</td></tr></table>	PRODUCT__ep2d_diff__p3_sms1
<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms1</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms1</td></tr></table>	PRODUCT__ep2d_diff__p3_sms1	
<table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms1</td></tr></table>	PRODUCT__ep2d_diff__p3_sms1		
PRODUCT__ep2d_diff__p3_sms1			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms2</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms2</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms2</td></tr></table>	PRODUCT__ep2d_diff__p3_sms2
<table border="1"><tr><td><table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms2</td></tr></table></td></tr></table>	<table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms2</td></tr></table>	PRODUCT__ep2d_diff__p3_sms2	
<table border="1"><tr><td>PRODUCT__ep2d_diff__p3_sms2</td></tr></table>	PRODUCT__ep2d_diff__p3_sms2		
PRODUCT__ep2d_diff__p3_sms2			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb1</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb1</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb1</td></tr></table>	C2P__cmrr_mbep2d_diff__p3_mb1
<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb1</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb1</td></tr></table>	C2P__cmrr_mbep2d_diff__p3_mb1	
<table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb1</td></tr></table>	C2P__cmrr_mbep2d_diff__p3_mb1		
C2P__cmrr_mbep2d_diff__p3_mb1			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb2</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb2</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb2</td></tr></table>	C2P__cmrr_mbep2d_diff__p3_mb2
<table border="1"><tr><td><table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb2</td></tr></table></td></tr></table>	<table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb2</td></tr></table>	C2P__cmrr_mbep2d_diff__p3_mb2	
<table border="1"><tr><td>C2P__cmrr_mbep2d_diff__p3_mb2</td></tr></table>	C2P__cmrr_mbep2d_diff__p3_mb2		
C2P__cmrr_mbep2d_diff__p3_mb2			

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\localizer

TA: 14 sec Coil Selection: Auto Voxel Size: 0.5×0.5×5.0 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	On
Auto Open Inline Display	On
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	All Segments
Inline Movie	Off

Routine

Slice Group	1
Slices	7
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	2
Slices	7
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	3
Slices	7
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	5.0 mm
TR	7.0 ms
TE	3.11 ms
Averages	1
Concatenations	13
AutoAlign	---

Contrast - Common

TR	7.0 ms
TE	3.11 ms
TD	0.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	10 deg

Contrast - Common

Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
SWI	Off
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	5.0 mm
Base Resolution	256
Phase Resolution	100 %
Interpolation	On

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	2
Reference Lines PE	24
Deep Resolve	Off
Phase Partial Fourier	Off
Asymmetric Echo	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	On
Distortion Correction	2D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slice Group	1
Slices	7
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	2
Slices	7
Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	3
Slices	7

Geometry - Common

Distance Factor	400 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	5.0 mm
TR	7.0 ms
Multi-Slice Mode	Sequential
Series	Interleaved
Concatenations	13

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	2
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slice Group	3
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
Special Saturation	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slice-sel.
LR Balancing	Off

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	7.0 ms
Segments	1
Concatenations	13

Physio - Cardiac

Tagging	None
Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	256 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	13

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1

Inline - Subtraction

StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Soft Tissue

Wash-in	Off
Wash-out	Off
TTP	Off
PEI	Off
MIP Time	Off
Measurements	1

Inline - Composing**Inline - MapIt**

MapIt	None
Flip Angle	10 deg
Measurements	1
Contrasts	1
TE	3.11 ms
TR	7.0 ms
Save Original Images	On

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	qfl
Dimension	2D
Excitation	Slice-sel.
RF Pulse Type	Normal
Gradient Mode	Whisper
Flow Compensation	None
Bandwidth	610 Hz/Px
Asymmetric Echo	Off
Segments	1

Sequence - Part 2

Introduction	On
RF Spoiling	On
Acoustic noise reduction	On

Sequence - Nuclei

TX/RX Nucleus	1H
---------------	----

Sequence - Nuclei

TX/RX Delta Frequency	0 Hz
TX Nucleus	None
TX Delta Frequency	0 Hz

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\PRODUCT_ep2d_bold_p3_sms1

TA: 11 sec Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	20.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	20.00 ms
MTC	Off
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	3
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	64

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	3
Reference Lines PE	30
SMS Factor	1
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On

BOLD

Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	3
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	2368 Hz/Px
Echo Spacing	0.55 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
--------------	-----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\PRODUCT_ep2d_bold_p3_sms2

TA: 9 sec Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc:: 6 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	20.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	20.00 ms
MTC	Off
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	3
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	64

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	3
Reference Lines PE	30
SMS Factor	2
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On

BOLD

Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	3
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	2368 Hz/Px
Echo Spacing	0.59 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
--------------	-----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\PRODUCT_ep2d_bold_p3_sms5

TA: 9 sec Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc:: 15 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	20.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	20.00 ms
MTC	Off
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	3
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	64

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	3
Reference Lines PE	30
SMS Factor	5
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On

BOLD

Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	3
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	2368 Hz/Px
Echo Spacing	0.69 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
--------------	-----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\C2P_cmrr_mbep2d_bold_p3_mb1

TA: 11 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	20.00 ms
Averages	1
Multi-band accel. factor	1
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	20.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	3
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	30
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Multi-band accel. factor	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0

BOLD

Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	3
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Dimension	2D
Excitation	Standard
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	2170 Hz/Px
Echo Spacing	0.59 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
RF Spoiling	Off

Sequence - Special

Excite pulse duration	2560 us
Min. prep scans	0
Delay before PC scans	0 us
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Sinc exc. pulse BWTP	5.20
Physio recording	Off
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\C2P_cmrr_mbep2d_bold_p3_mb2

TA: 22 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	20.00 ms
Averages	1
Multi-band accel. factor	2
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	20.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	3
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	30
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Multi-band accel. factor	2

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0

BOLD

Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	3
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Dimension	2D
Excitation	Standard
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	2170 Hz/Px
Echo Spacing	0.63 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
RF Spoiling	Off

Sequence - Special

Excite pulse duration	2560 us
Min. prep scans	0
Min. prep scans SB	0
Delay before PC scans	0 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
Opt. MB RF pulse BW	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Sinc exc. pulse BWTP	5.20
Physio recording	Off
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\C2P_cmrr_mbep2d_bold_p3_mb5

TA: 44 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
TE	20.00 ms
Averages	1
Multi-band accel. factor	5
AutoAlign	---

Contrast - Common

TR	1230.0 ms
TE	20.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	42 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	3
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	30
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1230.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	5

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	1230.0 ms
Multi-band accel. factor	5

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0

BOLD

Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	3
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Dimension	2D
Excitation	Standard
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	2170 Hz/Px
Echo Spacing	0.59 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
RF Spoiling	Off

Sequence - Special

Excite pulse duration	2560 us
Min. prep scans	0
Min. prep scans SB	0
Delay before PC scans	0 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
Opt. MB RF pulse BW	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Sinc exc. pulse BWTP	5.20
Physio recording	Off
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\PRODUCT__ep2d_diff__p3_sms1

TA: 36 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	80.00 ms
Concatenations	1
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	80.00 ms
MTC	Off
Magn. Preparation	None
Fat-Water Contrast	Fat Saturation
Fat Saturation	Weak
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	64

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	3
Reference Lines PE	30
SMS Factor	1
Deep Resolve	Off
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Concatenations	1

Physio - PACE

Resp. Control	Off
Concatenations	1

Diff

Diffusion Mode	MDDW
Diff. Directions	6
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	1000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	2368 Hz/Px
Echo Spacing	0.55 ms
Free Echo Spacing	Off
Optimization	None
EPI Factor	64

Sequence - Part 2

Introduction	Off
Phase Correction	Internal

Sequence - Assistant

SAR Assistant	Off
Optimization	None

\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\PRODUCT__ep2d_diff__p3_sms2

TA: 28 sec Coil Selection: Auto Voxel Size: 2.0×2.0×2.0 mm³ Acc:: 6 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	81.00 ms
Concatenations	1
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	81.00 ms
MTC	Off
Magn. Preparation	None
Fat-Water Contrast	Fat Saturation
Fat Saturation	Weak
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	64

Resolution - Common

Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	SMS
Reference Scans	EPI/Separate
Acceleration Factor PE	3
Reference Lines PE	30
SMS Factor	2
Deep Resolve	Off
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Concatenations	1

Physio - PACE

Resp. Control	Off
Concatenations	1

Diff

Diffusion Mode	MDDW
Diff. Directions	6
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	1000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Normal
Bandwidth	2368 Hz/Px
Echo Spacing	0.59 ms
Free Echo Spacing	Off
Optimization	None
EPI Factor	64

Sequence - Part 2

Introduction	Off
Phase Correction	Internal

Sequence - Assistant

SAR Assistant	Off
Optimization	None

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\C2P_cmrr_mbep2d_diff_p3_mb1

TA: 36 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	80.00 ms
Multi-band accel. factor	1
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	80.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	90 deg
Refocus flip angle	180 deg
Fat-Water Contrast	Fat Saturation
Grad. rev. fat suppr.	Enabled
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	30
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	1

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Rotation	0.00 deg
R >> L	128 mm
F >> H	128 mm
A >> P	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Multi-band accel. factor	1

Physio - PACE

Resp. Control	Off
Multi-band accel. factor	1

Diff

Diffusion Mode	MDDW
Diff. Directions	6
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	1000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Dimension	2D
Excitation	Standard
Gradient Mode	Normal
Bandwidth	2232 Hz/Px
Echo Spacing	0.57 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
RF Spoiling	Off
Phase Correction	Internal

Sequence - Special

Excite pulse duration	2560 us
Refocus pulse duration	5120 us
Min. prep scans	0
Delay before PC scans	0 us
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Physio recording	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\PHANTOM\PHANTOM\check_MBfactor_dcm2niix_issue870\C2P_cmrr_mbep2d_diff_p3_mb2

TA: 1:03 min Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
TE	80.00 ms
Multi-band accel. factor	2
AutoAlign	---

Contrast - Common

TR	3000.0 ms
TE	80.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle	90 deg
Refocus flip angle	180 deg
Fat-Water Contrast	Fat Saturation
Grad. rev. fat suppr.	Enabled
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Multiple Series	Off
Delay in TR	0.00 ms

Resolution - Common

FOV Read	128 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	2.0 mm
Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	30
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	10
Distance Factor	0 %
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
Phase Oversampling	0 %
FOV Read	128 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	3000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

Geometry - AutoAlign

Slice Group	1
Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P25.7 F12.9
L	0.0 mm
P	25.7 mm
F	12.9 mm
Initial Orientation	Coronal
Initial Rotation	0.00 deg

Geometry - Navigator

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 P25.7 F12.9 mm
Orientation	Coronal
Rotation	0.00 deg
R >> L	128 mm
F >> H	128 mm
A >> P	20 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.178714 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3000.0 ms
Multi-band accel. factor	2

Physio - PACE

Resp. Control	Off
Multi-band accel. factor	2

Diff

Diffusion Mode	MDDW
Diff. Directions	6
Diffusion Scheme	Bipolar
Diff. Weightings	2
b-value 1	0 s/mm ²
b-value 2	1000 s/mm ²
Averages 1	1
Averages 2	1
Dynamic Field Correction	Off
Invert Gray Scale	Off
Diff. Weighted Images	On
Trace Weighted Images	Off
Tensor	Off
FA Maps	Off
ADC Maps	Off
Exponential ADC Maps	Off
Calculated Image	Off

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Dimension	2D
Excitation	Standard
Gradient Mode	Normal
Bandwidth	2232 Hz/Px
Echo Spacing	0.61 ms
Free Echo Spacing	Off
EPI Factor	64

Sequence - Part 2

Introduction	Off
RF Spoiling	Off
Phase Correction	Internal

Sequence - Special

Excite pulse duration	2560 us
Refocus pulse duration	5120 us
Min. prep scans	0
Min. prep scans SB	0
Delay before PC scans	0 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
Opt. MB RF pulse BW	Off
Time-shifted MB RF	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off

Sequence - Special

Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Physio recording	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----